

# SYST'AM® P311C / GEL CUSHION

GEL CUSHION MADE OF VISCOELASTIC GEL

## MATERIALS



Visco Gel



Removable POLYMAILLE® cover



Removable POLYMAILLE® HD cover



### Cleaning of the gel:



(Do not immerse in water)

### Cleaning POLYMAILLE® cover:




### Cleaning POLYMAILLE® HD cover:



## INDICATIONS

### Risk





P311C / Gel cushion can be used, when:

- Seated position is necessary resulting from diseases or handicap and there is at the same time an increased risk of pressure sore.

### VISCOELASTIC HIGH-RETICULATED POLYURETHANE GEL

→ Avoids gel migration and prevents the bottoming effect of the bony protuberances.



100% IMPERMEABLE  
→ Completely adapted to incontinent patients.

## AVAILABLE VERSIONS





SYST'AM® P311C /Gel Cushion  
With POLYMAILLE® cover.



SYST'AM® P311C /Gel Cushion  
With POLYMAILLE® HD cover  
(with welded seams).



SYST'AM® P313C /Gel Cushion  
Gel-Cushion with a removable  
central pad.

| DESIGNATION                                     | ITEM CODE    | <br>(W x L - cm / inches) | <br>min - MAX<br>(kg / lbs) |
|---|--------------|---|--|
| SYST'AM® P311C<br>with POLYMAILLE® cover        | P311C40401HW | 40 x 40 cm / 16 x 16"   | 40 - 120 kg / 88 - 264 lbs   |
|   | P311C45431HW | 45 x 43 cm / 18 x 17"   | 40 - 140 kg / 88 - 308 lbs   |
| SYST'AM® P311C<br>with POLYMAILLE® HD cover     | P311C40401HF | 40 x 40 cm / 16 x 16"   | 40 - 120 kg / 88 - 264 lbs   |
|   | P311C45431HF | 45 x 43 cm / 18 x 17"   | 40 - 140 kg / 88 - 308 lbs   |
| SYST'AM® P313C with<br>a removable central pad. | P313C40401HW | 40 x 40 cm / 16 x 16"   | 40 - 100 kg / 88 - 220 lbs   |

