

ORDER FORM

BEN 4 plus





















creative in care

Release date 2024/11

Company	Invoice address	Order number / reference	
Street		Customer reference	
E-Mail		Date	
City, Zip code		Offer	
Telephone		Order	
Company	Delivery address		
Street	City, Zip code		

Standard equipment:	Size / order codes				
	MINI	STD	MAXI	SM 42	SM 45
<ul style="list-style-type: none"> • fully reclining backrest • seat with tilt-in space function (3 positions) • height variations for attaching of the fixation devices and accessories • angle and height adjustable footrest • 175-265 mm solid wheels • detachable grab bar with cover • height adjustable ergonomic handlebar 	front swivel castor-wheels /FSCW				
	B4P30RP00	B4P34RP00	B4P38RP00	B4P42RP00	B4P45RP00
	€	€	€	€	€

Seat width [mm]	200/step/300	240/step/340	280/step/380	320/step/420	350/step/450
Seat depth [mm]	220	290	370	370	460
Backrest height [mm]	540	640	710	710	840
Seat to footplate [mm]	140-280	200-330	230-400	230-400	240-460
Backrest angle	90° - 165°	90° - 165°	90° - 155°	90° - 155°	90° - 155°
Seat angle	19°/step/27°	20°/step/28°	19°/step/27°	18°/step/27°	20°/step/30°
Maximum load	40,0 kg	40,0 kg	50,0 kg	75,0 kg	75,0 kg

colour frame / pattern colour																									
frame	K1 - Anthracite																								
fabric patterns	<table border="0"> <tr> <td>R54</td> <td></td> <td>R57</td> <td></td> </tr> <tr> <td></td> <td>Slate grey</td> <td></td> <td>Antik pink</td> </tr> <tr> <td>R55</td> <td></td> <td>R58</td> <td></td> </tr> <tr> <td></td> <td>Denim blue</td> <td></td> <td>Papaya orange</td> </tr> <tr> <td>R56</td> <td></td> <td>R59</td> <td></td> </tr> <tr> <td></td> <td>Peas green</td> <td></td> <td>Ruby red</td> </tr> </table>	R54		R57			Slate grey		Antik pink	R55		R58			Denim blue		Papaya orange	R56		R59			Peas green		Ruby red
	R54		R57																						
		Slate grey		Antik pink																					
R55		R58																							
	Denim blue		Papaya orange																						
R56		R59																							
	Peas green		Ruby red																						

BEN 4 plus Accessories

Section E: Positioning and fixation - Head


No.	Desc. / Codes	B4P MINI	B4P STD	B4P MAXI	B4P SM 42	B4P SM 45	
E1	Headrest CLASSIC RPKE001	S	M	M	M	L	€
E2	Neckrest CLASSIC RPKE003	S	M	M	M	M	€
E3	Headwings BASIC RPKE007	S	M	M	M	L	€



Section G: Positioning and fixation - Trunk


G1	Side trunkrest CLASSIC RPKG001	S	M	M	L(XL)	XL	€
G2	Hip rest CLASSIC RPKG002	M	M	L	L	L	€
G3	Chest immo. belt CLASSIC RPKG013	S	M	M	L	L	€
G4	2-point immobilisation belt BASIC RPKG024	S	M	M	L	L	€
G5	4-point immobilisation belt BASIC RPKG025	S	M	M	L	L	€
G6	Check-belt addition to 4P-belt RPKG006	UNI	UNI	UNI	UNI	UNI	€
G7	Belt protector W25 RPKG007	UNI	UNI	UNI	UNI	UNI	€
G8	Belt protector W40 RPKG008	UNI	UNI	UNI	UNI	UNI	€
G9	5-point immobilisation belt BASIC RPKG026	S	M	M	L	L	€
G10	H-style harness BASIC RPKG027	S	M	M	L	L	€
G11	Backrest wedge - (reduction seat depth) RPKG015	S					€





Section I: Positioning and fixation - Pelvis / Abduction



No.	Desc. / Codes	B4P MINI	B4P STD	B4P MAXI	B4P SM 42	B4P SM 45	
I1	FLEXI-pads set (2H+1H+2H) RPKI003	S	M	M	M	M	€
I2	FLEXI-pads set with extended thigh guides RPKI004			UNI	UNI	UNI	€
I3	Abduction block BASIC RPKI012	S	M	M	L		€
I4	Abduction block CLASSIC RPKI005					L	€



Section K: Positioning and fixation - Arms


K1	Padded frame cover RPKK001	UNI	UNI	UNI	UNI	UNI	€
K2	Extended side protection CLASSIC RPKK002	UNI	UNI	UNI	UNI		€
K3	Padded frame cover with s.p. BASIC RPKK004	M	M	M	L	L	€
K4	Grab bar RNDK005	UNI(SE)	UNI(SE)	UNI(SE)	UNI(SE)	UNI(SE)	€
K5	Armrest SOFT RPKK005	S	S	M	M	L	€
K6	Armrest CLASSIC adj. RPKK00630	UNI					€
K7	Armrest CLASSIC adj. RPKK00634		UNI				€
K8	Armrest CLASSIC adj. RPKK00638			UNI			€
K9	Armrest CLASSIC adj. RPKK00642				UNI		€

K10	Armrest CLASSIC adj. RPKK00645					UNI	€
K11	Work board BASIC RPKK111	UNI	UNI	UNI			€
K12	Work board BASIC RPKK112				UNI	UNI	€







Section M: Positioning and fixation - Legs

M1	Feet fixation - pair RPKM001	UNI	UNI	UNI	UNI	UNI	€
M2	Footplate soft padding RPKM002	M	M	L	L	L	€
M3	Calf rest (textile) RPKM004	UNI	UNI	UNI	UNI	UNI	€



Section P: Protection and comfort

P1	Canopy BASIC RPKP003	S	M	M	L	L	€
P2	Canopy CLASSIC RPKP001	(S)	(M)	(M)	(L)	(L)	€
P3	Rain cover RPKP051	M	M	L	L	XL	€
P4	Footmuff - Summer RPKP101	S	M	M	L	L	€
P5	Footmuff - Winter RPKP151	S	M	M	L	L	€






Section R: Bags, basket

No.	Desc. / Codes	B4P MINI	B4P STD	B4P MAXI	B4P SM 42	B4P SM 45	
R1	Bag JUMBO RPKR001	UNI	UNI	UNI	UNI	UNI	€
R2	Backpack FOLD RPKR002	UNI	UNI	UNI	UNI	UNI	€
R2	Basket (textile) RPKR051	M	M	L	XL	XL	€
R4	Basket-HD (textile) RPKR101	M	M	L	XL	XL	€










Section T: Accessory - specialities

T1	Backrest HD-telescope w.tilt-lock RPKT004	UNI	UNI				€
T2	Backrest HD-telescope w.tilt-lock RPKT005			UNI	UNI	UNI	€



Section U: Accessory - set

No.	Desc. / Codes	B4P MINI	B4P STD	B4P MAXI	B4P SM 42	B4P SM 45	
U1	Air Rear tyre (set) (RW only) RPKU001	UNI	UNI	UNI	UNI	UNI	€
U2	Air Rear tyre (set) (FFW only) RPKU002	(UNI)	(UNI)	(UNI)			€
U3	Front Fork with swivel lock (set) RPKU051	UNI	UNI	UNI	UNI	UNI	€
U4	Footrest-SHORT version RPKU20101	UNI					€
U5	Footplate angle adjustable RPKU20201	UNI	UNI				€
U6	Footplate angle adjustable RPKU20202			UNI	UNI	UNI	€



Sizes:

XS - extra small
 S - small
 M - medium
 L - large
 XL - extra large
 UNI - universal size (one size only)

Other notes:

(SE) Standard Equipment
 (empty position) - accessory is not available