

# ORDER FORM

# TAMPA





**patron**<sup>®</sup>  
creative in care

Release date 2023/12

Company	Invoice address	Order number / reference	
Street		Customer reference	
E-Mail		Date	
City, Zip code		Offer	
Telephone		Order	
Company	Delivery address		
Street	City, Zip code		

Standard equipment:	MINI	STD	MAXI	SM 42
<ul style="list-style-type: none"> <li>• seat depth adjustment</li> <li>• height adjustable backrest</li> <li>• width adjustable seat by FlexiPads set (optimal equipment)</li> <li>• seat height and tilt in space function</li> <li>• angle and height adjustable footrest</li> <li>• height adjustable armrest</li> <li>• rear castor wheels with brake</li> <li>• telescopic wheeled base with length adjustment</li> <li>• grab bar, armrest, push handle.</li> </ul>	front swivel castor-wheels /FSCW			
	TMC30AM00	TMC34AM00	TMC38AM00	TMC42AM00
	€	€	€	€

Backrest height [mm]	450-600	520-680	580-680	580-690
Seat depth [mm]	200-270	280-310	320-430	320-430
Seat width [mm]	200/step/300	240/step/340	280/step/380	320/step/420
Seat to footrest [mm]	140-290	190-340	240-390	240-390
Backrest angle	90° - 165°	90° - 165°	90° - 160°	90° - 160°
Seat angle	+/-27°	+/-27°	+/-25°	+/-25°
Maximum load	40,0 kg	40,0 kg	50,0 kg	75,0 kg

colour frame / pattern colour	
<b>frame</b>	<b>K1 - Anthracite</b>
<b>fabric patterns</b>	R54  Slate grey <span style="margin-left: 200px;">R57  Antik pink</span>
	R55  Denim blue
	R56  Peas green

## TAMPA Accessories

### Section E: Positioning and fixation - Head


No.	Desc. / Codes	TMC MINI	TMC STD	TMC MAXI	TMC SM42	€
E1	Headrest CLASSIC RPKE001	S	M	M	M	€
E2	Headwings CLASSIC RPKE002	(XS ) S	M	M	M	€
E3	Neckrest CLASSIC RPKE003	S	M	M	M	€
E4	Headwings RPKE004				L	€
E5	Backrest ext.for Headwings RPKE005				UNI	€
E6	Headrest contoured, height+angle adj. RPKE006				L	€



### Section G: Positioning and fixation - Trunk


No.	Desc. / Codes	TMC MINI	TMC STD	TMC MAXI	TMC SM42	€
G1	Side trunkrest CLASSIC RPKG001	S	M	M	(L )XL	€
G2	Hip rest CLASSIC RPKG002	M	M	L	L	€
G3	Chest immo. belt CLASSIC RPKG0013	S	M	M	L	€
G4	2-point immobilisation belt CLASSIC RPKG003	S	M	M	M	€
G5	4-point immobilisation belt CLASSIC RPKG004	S	M	M	L	€
G6	Chest-belt addition to 4P-belt RPKG006	UNI	UNI	UNI	UNI	€
G7	Belt protector W25 (belt 25mm) RPKG007	UNI	UNI	UNI	UNI	€
G8	Belt protector W40 (belt 40mm) RPKG008	UNI	UNI	UNI	UNI	€
G9	5-point immobilisation belt CLASSIC RPKG005	S	M	M	L	€
G10	Immo. Waistcoat CLASSIC RPKG009	(XS ) S	M	M	M	€
G11	H-style harness CLASSIC RPKG011	(XS ) S	M	M	L	€
G12	Immo. Waistcoat TBX RPKG010	(XS ) S	M	M	M	€

G13 H-style harness TBX RPKG012	(XS ) S	M	M	M	€
G14 Lumbar support W15, 25, 35 RPKG014	S	M	L	L	€
G15 Backrest wedge RPKG015	S				€



### Section I: Positioning and fixation - Pelvis / Abduction

No.	Desc. / Codes	TMC MINI	TMC STD	TMC MAXI	TMC SM42	
I1	Pelvis fixation RPKI001	S	S	M	M	€
I2	Abduction belt CLASSIC RPKI002	(XS ) S	M	M	M	€
I3	FLEXI-pads set RPKI003	S	M	M	M	€
I4	FLEXI-pads set with extended thigh guides RPKI004			UNI	UNI	€
I5	Abduction block CLASSIC RPKI005	XS	S	M	M	€





### Section K: Positioning and fixation - Arms

No.	Desc. / Codes	TMC MINI	TMC STD	TMC MAXI	TMC SM42	
K1	Grab bar	UNI(SE)	UNI(SE)	UNI(SE)	UNI(SE)	€
K2	Work board CLASSIC RPKK101	UNI	UNI	UNI	UNI	€
K3	Work board BASIC RPKK111	UNI	UNI	UNI		€
K4	Work board BASIC - XL RPKK112				UNI	€





### Section M: Positioning and fixation - Legs

No.	Desc. / Codes	TMC MINI	TMC STD	TMC MAXI	TMC SM42	
M1	Feet fixation - pair RPKM001	UNI	UNI	UNI	UNI	€
M2	Footplate soft padding RPKM002	M	M	L	L	€
M3	Calf rest (textile) RPKM004	UNI	UNI	UNI	UNI	€



## Section U: Accessory - set

No.	Desc. / Codes	TMC MINI	TMC STD	TMC MAXI	TMC SM42	
U1	Footrest - SHORT version RPKU20101	UNI				€
U2	Footplate angle adjustable RPKU20201	UNI	UNI			€
U3	Footplate angle adjustable RPKU20202			UNI	UNI	€



### Sizes:

XS - extra small  
 S - small  
 M - medium  
 L - large  
 XL - extra large  
 UNI - universal size (one size only)

### Other notes:

( SE ) Standard Equipment  
 (empty position) - accessory is not available